












## A-TYPE MICRO MATIC KEG COUPLER SERVICE TIPS SERVICE KIT # 35776711

 <p>1. Unclip handle hinge pin (#35654034) &amp; slide off. Slide handle (#35654022) back to remove.</p>	 <p>2. Push Probe down partially and remove main base seal (large black one - 35700003)</p>	 <p>3. Push the probe through the body and remove from the bottom of coupler.</p>
 <p>4. There are 2x black Seals (#35700004) in the body of the keg coupler. Remove these and replace. Use small amount Kluber grease on O rings</p>		 <p>5. Replace the probe by inserting the threaded end of the probe into the body and slide in as shown below. Refit new keg seal (#35700003)</p>
	 <p>7. Line the probe 'flats' in line with hinge pin entry. Push in hinge pin and fasten in place.</p>	 <p><b>A Type Service Kit - 35776711</b> <b>Kluber Grease - 47517171</b> <b>O-Ring Pick Set - 00000090</b></p>